

Identification of Fungicide –Thiram using SERS

INTRODUCTION

- > Fungicides are used for high productivity in farming
- > Retained fungicides in soil/water is harmful
- > Existing detection methods are time consuming and sophisticated
- > Surface-enhanced Raman spectroscopy (SERS) – high sensitivity, short detection time and nondestructive detection
- > In the current study, Thiram has been detected using SERS with 10^3 enhancement


EXPERIMENTAL CONDITION

- > Raman Spectrometer:
- > Model: IndiRAM P785_Tracer
- > Laser excitation source
 - Wavelength: 785 nm
 - Power: 360 mW
- > Exposure time: 5 sec
- > No of Acquisition: 5
- > SERS: Liquid SERS


CONCLUSION

The vibration modes of Thiram were quickly identified by using SERS Technique.

RESULT & DISCUSSION



SERS



- >SERS shown enhanced Raman Signal
- >Detected minimum of 0.1 μ M of Thiram on varying Concentration


REFERENCES

1. T.C. Dao, T.Q.N. Luong, T.A. Cao, N.M. Kieu, Advances in Natural Sciences: Nanoscience and Nanotechnology, 10 (2019) 025012.
2. H. Pu, Z. Huang, F. Xu, D.-W. Sun, Food Chemistry, 343 (2021) 128548.

OUR PRODUCTS



IndiRAM™ CTR Series



IndiRAM™ CTR for
Quantum Characterizer




IndiRAM™ CTR Series
Chiral Raman Spectrometer



IndiRAM™ Edu



IndiRAM™ CTR
P-Series



IndiRAM™ Portable
System



Diamond Detection
System



IndiRAM™
Handy Series



IndiRAM™ CTR-
Mini Series



Solutions through Technology



G1-43, RIICO Industrial Area, Sitapura, Jaipur, Rajasthan
302022



www.technosphotonics.com



info@technos.in



+91 8003220052